



Comfort, Health and Style!

PRODUCT CATALOGUE

BUILDING A HEALTHY COMMUNITY
AT HOME

THERAPIES. HAND IN HAND.



JOBST® INNOVATIONS

BSN medical connects its broad portfolio in Wound / Vascular Management, Lymphology and non-invasive Orthopedics with strong technological competencies to meet the market's need for reduced complexity and more efficient treatments. With a continuously developing unique and holistic network of integrated therapy solutions, BSN medical consistently supports a continuum of care.



JOBST® SOFTFIT™ – True Gradient Compression

The SoftFit™ Technology is a unique silicone yarn that is knitted into the inner knee welt of the stocking. It adheres softly to the skin, increases the gripping properties and keeps the stockings in place. This innovative grip zone reduces constriction and therefore offers **true gradient compression**.



3D KNITTING – Like a Second Skin

The 3D knitting technology reproduces each body shape accurately and has revolutionised the market with compression garments of uncompromising quality and an anatomical fit that feels **like a second skin**.



JOBST® SENSITIVE – More Comfortable Silicone

JOBST® has enhanced its high quality thigh high garments with a **more comfortable silicone** band with yarns that are woven into the welt fabric, not coated onto the material. This allows for even distribution of friction, soft grip on the skin, less trapped moisture which leads to a reduced risk for red marks and irritated skin.



JOBST® ADVANCED COMFORT – Superior Wearing Comfort

The new JOBST® Advanced Comfort system consists of the unique micro fiber yarns to provide patients with **superior wearing comfort**. The micro fiber yarn construction and the finishes enable the fabric to take up moisture easily and spread it fast over a wide surface area.



DRI-RELEASE® – Dry, Cool and Comfortable Feet

Superior moisture management to help keep feet **dry, cool and comfortable**.



C² EFFECT – Natural Odor and Moisture Control

The unique combination of cotton & carbon provide high wearing comfort. The premium quality cotton ensures a natural feeling to the skin thus increasing comfort. Carbon, in its activated form, supports **natural odor and moisture control**.

TABLE OF CONTENTS

JOBST® Innovations

3

VASCULAR

READY-TO-WEAR STOCKINGS

Styles and Colours Overview	6
Sizing Charts	8
JOBST® UltraSheer	9
JOBST® Opaque	16
JOBST® soSoft	20
JOBST® Sport	21
JOBST® ActiveWear	22
JOBST® Travel Socks	23
JOBST® Relief	24
JOBST® forMen Ambition	26
JOBST® forMen Casual	29
JOBST® forMen	30

CUSTOM STOCKINGS

JOBST® Custom Seamed	32
JOBST® Custom Seamless Soft	36
JOBST® Bellavar®	38

LYMPHOLOGY

CUSTOM GARMENTS

JOBST® Elvarex®	40
JOBST® Elvarex® Soft	41
JOBST® Elvarex® Plus	42
JOBST® Relax	43
Lymphology Custom Garments Options	44

READY-TO-WEAR GARMENTS

JOBST® Bella™ Lite	45
JOBST® Bella™ Strong	46

TABLE OF CONTENTS

BANDAGES/WRAPS

JOBST® FarrowWrap®	48
Comprilan®	54
Comprifoam®	55
Artiflex®	55
JOBST® Compri2	56
JOBST® UlcerCARE™	57

OTHERS

JOBST® Surgical vest	58
Stockings Accessories & Donning Aids	59
Fiber Contents	60

All products in this catalogue are Latex Free, except : JOBST® Elvarex®, Comprilan® (clips only), Foam rubber pad and donning gloves.



STYLES AND COLOURS OVERVIEW

WOMEN PRODUCTS

Style	Knee high closed toe	Knee high open toe	Knee high SoftFit™ closed toe	Thigh high closed toe	Thigh high open toe	Pantyhose	Maternity	Petite
15-20 mmHg								
JOBST® UltraSheer								
JOBST® Opaque								
JOBST® soSoft								
20-30 mmHg								
JOBST® UltraSheer								
JOBST® Opaque								
JOBST® soSoft								
30-40 mmHg								
JOBST® UltraSheer								
JOBST® Opaque								
JOBST® soSoft								

STYLES AND COLOURS OVERVIEW

UNISEX PRODUCTS

Style	Knee high closed toe	Knee high open toe	Thigh high closed toe	Thigh high open toe	Pantyhose
15-20 mmHg					
JOBST® Sport					
JOBST® ActiveWear					
JOBST® Travel Socks					
JOBST® Relief					
20-30 mmHg					
JOBST® Sport					
JOBST® ActiveWear					
JOBST® Relief					
30-40 mmHg					
JOBST® ActiveWear					
JOBST® Relief					

MEN PRODUCTS



Style	Knee high closed toe	Knee high open toe	Knee high SoftFit™ closed toe	Thigh high open toe	Long
15-20 mmHg					
JOBST® forMen Ambition	  		  		  
JOBST® forMen Casual	 				
JOBST® forMen	  			 	
20-30 mmHg					
JOBST® forMen Ambition	  		  		  
JOBST® forMen Casual	 				
JOBST® forMen	 			 	
30-40 mmHg					
JOBST® forMen Ambition	  		  		  
JOBST® forMen Casual	 				
JOBST® forMen	 			 	

SIZING CHARTS

JOBST® (all except products below)

Measurements (cm)	S	M	L	XL	LFC*	XLFC**
Ankle (B)	18-21	21-25	25-29	29-33	25-30	30-36
Calf (C)	28-38	30-42	32-46	34-50	46-61	46-61
Thigh (G)	40-62	46-70	54-78	60-81		
Hip (H)	71-117	76-127	81-137	102-166		
Hip (Maternity)	81-132	86-142	91-152	112-180		
Women Knee high - Regular	≥ 38					
Women Knee high - Petite	< 38					
Women Thigh high - Regular	≥ 68.5					
Women Thigh high - Petite	< 68.5					
Men - Regular	≤ 50					
Men - Tall	> 50					

Measuring Points



JOBST® UltraSheer Aurora

Measurements (cm)	1	2	3	4	5	6
Ankle (B)	18-20	20-22	22-24	24-26	26-28	28-30
Ankle and calf mid-point (B1)	23-27	24-29	26-32	29-35	31-37	33-39
Largest Calf (C)	28-34	30-37	33-40	35-43	37-46	39-49
Knee (D)	27-33	29-36	32-39	34-42	36-45	38-48
Thigh (G)	43-48	45-52	49-56	53-60	56-64	60-68
Waist (T)	80-89	90-99	90-99	110-119	110-119	120 +
Hips (H)	95-104	105-124	105-124	125-134	125-134	135 +
Women Knee high - Regular (ID)	39-45					
Women Thigh high - Regular (IG)	70-80					
Women Pantyhose - Regular (IK)	75-85					

JOBST® forMen Ambition

Measurements (cm)	1	2	3	4	5	6
Ankle (B)	18-20	20-22	22-24	24-27	27-30	30-34
Calf (C)	28-36	30-40	34-44	38-50	44-56	50-62
Regular	38-44					
Long	45-51					

JOBST® Travel Socks

Measurements (cm)	1	2	3	4	5
Ankle (B)	18-22	20-23	22-25	24-27	25-29
Calf (C)	27-37	31-41	33-43	36-46	40-50
Shoe size women	4.5 - 6.5	7 - 9	9.5 - 11	11.5 - 13	
Shoe size men	3.5 - 5.5	6 - 8	8.5 - 10	10.5 - 12	12.5 - 14

*LFC : Large Full Calf **LFC : Extra-Large Full Calf

JOBST® ULTRASHEER

SHEER ELEGANCE

- Unparalleled sheerness in a lightweight design
- The pattern design and trendy colours combine therapeutic efficacy and fashionable appearance


COMFORT & DURABILITY

- Reinforced heel and toe for greater durability
- High stretch capabilities for easier donning and better comfort

STYLES




TOE VERSIONS



COLOURS



BAND VARIATIONS



COMPRESSION CLASSES

15-20 mmHg

20-30 mmHg

30-40 mmHg

Not all colours and styles may be available in all compression ranges; please refer to product charts.

JOBST® UltraSheer Knee High



Colour	Style/Band	Toe Style	S	M	L	XL	LFC	XLFC
15-20 mmHg								
Black		Closed	7516504	7516604	7516704	7516804	7541601	7541701
	Pattern	Closed	7363300	7363301	7363302	7363303		
	Petite	Closed	7540801	7540901	7541001	7541501		
	SoftFit™	Closed	7768206	7768207	7768208	7768209	7768210	7768211
		Open	7516507	7516607	7516707	7516807		
Natural		Closed	7516525	7516625	7516725	7516825	7541602	7541702
	Petite	Closed	7540802	7540902	7541002	7541502		
	SoftFit™	Closed	7768200	7768201	7768202	7768203	7768204	7768205
		Open	7516530	7516630	7516730	7516830		
Honey		Closed	7616528	7516628	7516728	7516828		
Suntan		Closed	7516526	7516626	7516726	7516826		
	Petite	Closed	7540803	7540903	7541003	7541503		
Espresso		Closed	7516529	7516629	7516729	7516829		
	Pattern	Closed	7279600	7279601	7279602	7279603		
Anthracite		Closed	7516527	7516627	7516727	7516827		
Sun Bronze		Closed	7516531	7516631	7516731	7516831		
20-30 mmHg								
Black		Closed	7520104	7520204	7520304	7520404	7543301	7543401
	Pattern	Closed	7363400	7363401	7363402	7363403		
	Petite	Closed	7543501	7543601	7543701	7543801		
	SoftFit™	Closed	7768218	7768219	7768220	7768221	7768222	7768223
		Open	7520127	7520227	7520327	7520427	7543327	7543427
	Petite	Open	7543504	7543604	7543704	7543804		
Natural		Closed	7520119	7520219	7520319	7520419	7543302	7543402
	Petite	Closed	7543502	7543602	7543702	7543802		
	SoftFit™	Closed	7768212	7768213	7768214	7768215	7768216	7768217
		Open	7520129	7520229	7520329	7520429	7543329	7543429
	Petite	Open	7543505	7543605	7543705	7543805		
Honey		Closed	7520123	7520223	7520323	7520423		
		Open	7520128	7520228	7520328	7520428	7543328	7543428
Suntan		Closed	7520120	7520220	7520320	7520420		
	Petite	Closed	7543503	7543603	7543703	7543803		
Espresso		Closed	7520121	7520221	7520321	7520421		
	Pattern	Closed	7279700	7279701	7279702	7279703		
Anthracite		Closed	7520122	7520222	7520322	7520422		
Sun Bronze		Closed	7520130	7520230	7520330	7520430		
30-40 mmHg								
Black		Closed	7524104	7524204	7524304	7524404	7544701	7544801
	Petite	Closed	7544901	7545001	7545101	7545201		
	SoftFit™	Closed	7768230	7768231	7768232	7768233	7768234	7768235
		Open	7524127	7524227	7524327	7524427		
Natural		Closed	7524121	7524221	7524321	7524421	7544702	7544802
	Petite	Closed	7544902	7545002	7545102	7545202		
	SoftFit™	Closed	7768224	7768225	7768226	7768227	7768228	7768229
		Open	7524128	7524228	7524328	7524428		
Honey		Closed	7524124	7524224	7524324	7524424		
Suntan		Closed	7524122	7524222	7524322	7524422		
	Petite	Closed	7544903	7545003	7545103	7545203		
Espresso		Closed	7524125	7524225	7524325	7524425		
Anthracite		Closed	7524123	7524223	7524323	7524423		
Sun Bronze		Closed	7524129	7524229	7524329	7524429		

JOBST® UltraSheer Thigh High



Colour	Style/Band	Toe Style	S	M	L	XL
15-20 mmHg						
● Black	Lace	Closed	7516904	7517004	7517104	7517204
	Dot	Closed	7542501	7542601	7542701	7542801
	Pattern/Dot	Closed	7363500	7363501	7363502	7363503
	Petite/Lace	Closed	7542901	7543001	7543101	7543201
	Sensitive	Closed	7768604	7768605	7768606	7768607
	Petite/Sensitive	Closed	7768620	7768621	7768622	7768623
	Dot	Open	7542504	7542604	7542704	7542804
● Natural	Lace	Closed	7516922	7517022	7517122	7517222
	Dot	Closed	7542502	7542602	7542702	7542802
	Petite/Lace	Closed	7542902	7543002	7543102	7543202
	Sensitive	Closed	7768612	7768613	7768614	7768615
	Petite/Sensitive	Closed	7768628	7768629	7768630	7768631
	Dot	Open	7542505	7542605	7542705	7542805
● Honey	Lace	Closed	7516925	7517025	7517125	7517225
● Suntan	Lace	Closed	7516923	7517023	7517123	7517223
	Dot	Closed	7542503	7542603	7542703	7542803
	Petite/Lace	Closed	7542903	7543003	7543103	7543203
● Espresso	Lace	Closed	7516926	7517026	7517126	7517226
	Pattern/Dot	Closed	7279604	7279605	7279606	7279607
● Anthracite	Lace	Closed	7516924	7517024	7517124	7517224
● Sun Bronze	Lace	Closed	7516927	7517027	7517127	7517227
	Dot	Closed	7542506	7542606	7542706	7542806
20-30 mmHg						
● Black	Lace	Closed	7520504	7520604	7520704	7520804
	Dot	Closed	7543901	7544001	7544101	7544201
	Pattern/Dot	Closed	7363600	7363601	7363602	7363603
	Petite/Lace	Closed	7544301	7544401	7544501	7544601
	Sensitive	Closed	7768704	7768705	7768706	7768707
	Petite/Sensitive	Closed	7768720	7768721	7768722	7768723
	Dot	Open	7543904	7544004	7544104	7544204
● Natural	Lace	Closed	7520523	7520623	7520723	7520823
	Dot	Closed	7543902	7544002	7544102	7544202
	Petite/Lace	Closed	7544302	7544402	7544502	7544602
	Sensitive	Closed	7768712	7768713	7768714	7768715
	Petite/Sensitive	Closed	7768728	7768729	7768730	7768731
	Dot	Open	7543906	7544006	7544106	7544206
	Petite/Dot	Open	7543907	7544007	7544107	7544207
● Honey	Lace	Closed	7520527	7520627	7520727	7520827
	Dot	Open	7543905	7544005	7544105	7544205
● Suntan	Lace	Closed	7520524	7520624	7520724	7520824
	Dot	Closed	7543903	7544003	7544103	7544203
	Petite/Lace	Closed	7544303	7544403	7544503	7544603
● Espresso	Lace	Closed	7520525	7520625	7520725	7520825
	Pattern/Dot	Closed	7279704	7279705	7279706	7279707
● Anthracite	Lace	Closed	7520526	7520626	7520726	7520826
● Sun Bronze	Lace	Closed	7520531	7520631	7520731	7520831
	Dot	Closed	7543908	7544008	7544108	7544208


JOBST® UltraSheer Thigh High



Colour	Style/Band	Toe Style	S	M	L	XL
30-40 mmHg						
● Black	Lace	Closed	7524504	7524604	7524704	7524804
	Dot	Closed	7545301	7545401	7545501	7545601
	Petite/Lace	Closed	7545701	7545801	7545901	7546001
	Sensitive	Closed	7768804	7768805	7768806	7768807
	Petite/Sensitive	Closed	7768820	7768821	7768822	7768823
● Natural	Lace	Closed	7524520	7524620	7524720	7524820
	Dot	Closed	7545302	7545402	7545502	7545602
	Petite/Lace	Closed	7545702	7545802	7545902	7546002
	Sensitive	Closed	7768812	7768813	7768814	7768815
	Petite/Sensitive	Closed	7768828	7768829	7768830	7768831
● Honey	Dot	Open	7545304	7545404	7545504	7545604
	Lace	Closed	7524523	7524623	7524723	7524823
● Suntan	Lace	Closed	7524521	7524621	7524721	7524821
	Dot	Closed	7545303	7545403	7545503	7545603
	Petite/Lace	Closed	7545703	7545803	7545903	7546003
● Espresso	Lace	Closed	7524524	7524624	7524724	7524824
● Anthracite	Lace	Closed	7524522	7524622	7524722	7524822
● Sun Bronze	Lace	Closed	7524525	7524625	7524725	7524825
	Dot	Closed	7545305	7545405	7545505	7545605



JOBST® UltraSheer Pantyhose



Colour	Style/Band	Toe Style	S	M	L	XL
15-20 mmHg						
● Black		Closed	7517304	7517404	7517504	7517604
	Pattern	Closed	7363700	7363701	7363702	7363703
	Maternity	Closed	7517701	7517801	7517901	7518001
● Natural		Closed	7517312	7517412	7517512	7517612
	Maternity	Closed	7517704	7517804	7517904	7518004
● Honey		Closed	7517315	7517415	7517515	7517615
● Suntan		Closed	7517313	7517413	7517513	7517613
● Espresso		Closed	7517316	7517416	7517516	7517616
	Pattern	Closed	7279608	7279609	7279610	7279611
● Anthracite		Closed	7517314	7517414	7517514	7517614
● Sun Bronze		Closed	7517317	7517417	7517517	7517617
20-30 mmHg						
● Black		Closed	7520904	7521004	7521104	7521204
	Pattern	Closed	7363800	7363801	7363802	7363803
	Maternity	Closed	7521303	7521403	7521503	7521603
● Natural		Closed	7520913	7521013	7521113	7521213
	Maternity	Closed	7521306	7521406	7521506	7521606
● Honey		Closed	7520917	7521017	7521117	7521217
● Suntan		Closed	7520914	7521014	7521114	7521214
● Espresso		Closed	7520915	7521015	7521115	7521215
	Pattern	Closed	7279708	7279709	7279710	7279711
● Anthracite		Closed	7520916	7521016	7521116	7521216
● Sun Bronze		Closed	7520918	7521018	7521118	7521218
30-40 mmHg						
● Black		Closed	7524904	7525004	7525104	7525204
● Natural		Closed	7524913	7525013	7525113	7525213
● Honey		Closed	7524916	7525016	7525116	7525216
● Suntan		Closed	7524914	7525014	7525114	7525214
● Espresso		Closed	7524917	7525017	7525117	7525217
● Anthracite		Closed	7524915	7525015	7525115	7525215
● Sun Bronze		Closed	7524918	7525018	7525118	7525218

JOBST® UltraSheer Aurora Knee High



Colour	1	2	3	4	5	6
Class 1 (18-21 mmHg)						
Cranberry	73575-00048-00	73575-00049-00	73575-00050-00	73575-00051-00	73575-00052-00	73575-00053-00
Navy Blue	73575-00054-00	73575-00055-00	73575-00056-00	73575-00057-00	73575-00058-00	73575-00059-00
Lavender	73575-00060-00	73575-00061-00	73575-00062-00	73575-00063-00	73575-00064-00	73575-00065-00
Taupe	73575-00090-00	73575-00091-00	73575-00092-00	73575-00093-00	73575-00094-00	73575-00095-00
Crystal Grey	73575-00096-00	73575-00097-00	73575-00098-00	73575-00099-00	73575-00100-00	73575-00101-00
Cherry Red	73575-00102-00	73575-00103-00	73575-00104-00	73575-00105-00	73575-00106-00	73575-00107-00
Class 2 (22-33 mmHg)						
Cranberry	73583-00066-00	73583-00067-00	73583-00068-00	73583-00069-00	73583-00070-00	73583-00071-00
Navy Blue	73583-00072-00	73583-00073-00	73583-00074-00	73583-00075-00	73583-00076-00	73583-00077-00
Lavender	73583-00078-00	73583-00079-00	73583-00080-00	73583-00081-00	73583-00082-00	73583-00083-00
Taupe	73583-00126-00	73583-00127-00	73583-00128-00	73583-00129-00	73583-00130-00	73583-00131-00
Crystal Grey	73583-00132-00	73583-00133-00	73583-00134-00	73583-00135-00	73583-00136-00	73583-00137-00
Cherry Red	73583-00138-00	73583-00139-00	73583-00140-00	73583-00141-00	73583-00142-00	73583-00143-00

JOBST® UltraSheer Aurora Thigh High



Colour	1	2	3	4	5	6
Class 1 (18-21 mmHg)						
Cranberry	73577-00774-00	73577-00775-00	73577-00776-00	73577-00777-00	73577-00778-00	73577-00779-00
Navy Blue	73577-00714-00	73577-00715-00	73577-00716-00	73577-00717-00	73577-00718-00	73577-00719-00
Lavender	73577-00786-00	73577-00787-00	73577-00788-00	73577-00789-00	73577-00790-00	73577-00791-00
Taupe	73577-00738-00	73577-00739-00	73577-00740-00	73577-00741-00	73577-00742-00	73577-00743-00
Crystal Grey	73577-00798-00	73577-00799-00	73577-00800-00	73577-00801-00	73577-00802-00	73577-00803-00
Cherry Red	73577-00810-00	73577-00811-00	73577-00812-00	73577-00813-00	73577-00814-00	73577-00815-00
Class 2 (22-33 mmHg)						
Cranberry	73585-00802-00	73585-00803-00	73585-00804-00	73585-00805-00	73585-00806-00	73585-00807-00
Navy Blue	73585-00742-00	73585-00743-00	73585-00744-00	73585-00745-00	73585-00746-00	73585-00747-00
Lavender	73585-00814-00	73585-00815-00	73585-00816-00	73585-00817-00	73585-00818-00	73585-00819-00
Taupe	73585-00766-00	73585-00767-00	73585-00768-00	73585-00769-00	73585-00770-00	73585-00771-00
Crystal Grey	73585-00826-00	73585-00827-00	73585-00828-00	73585-00829-00	73585-00830-00	73585-00831-00
Cherry Red	73585-00838-00	73585-00839-00	73585-00840-00	73585-00841-00	73585-00842-00	73585-00843-00

Note : Delivery delays 5 to 7 business days for JOBST® UltraSheer® Aurora Edition.

JOBST® UltraSheer Aurora Pantyhose



Colour	1	2	3	4	5	6
Class 1 (18-21 mmHg)						
Cranberry	73581-00120-00	73581-00121-00	73581-00122-00	73581-00123-00	73581-00124-00	73581-00125-00
Navy Blue	73581-00126-00	73581-00127-00	73581-00128-00	73581-00129-00	73581-00130-00	73581-00131-00
Lavender	73581-00156-00	73581-00157-00	73581-00158-00	73581-00159-00	73581-00160-00	73581-00161-00
Taupe	73581-00282-00	73581-00283-00	73581-00284-00	73581-00285-00	73581-00286-00	73581-00287-00
Crystal Grey	73581-00294-00	73581-00295-00	73581-00296-00	73581-00297-00	73581-00298-00	73581-00299-00
Cherry Red	73581-00306-00	73581-00307-00	73581-00308-00	73581-00309-00	73581-00310-00	73581-00311-00
Class 2 (22-33 mmHg)						
Cranberry	73589-00174-00	73589-00175-00	73589-00176-00	73589-00177-00	73589-00178-00	73589-00179-00
Navy Blue	73589-00180-00	73589-00181-00	73589-00182-00	73589-00183-00	73589-00184-00	73589-00185-00
Lavender	73589-00246-00	73589-00247-00	73589-00248-00	73589-00249-00	73589-00250-00	73589-00251-00
Taupe	73589-00498-00	73589-00499-00	73589-00500-00	73589-00501-00	73589-00502-00	73589-00503-00
Crystal Grey	73589-00516-00	73589-00517-00	73589-00518-00	73589-00519-00	73589-00520-00	73589-00521-00
Cherry Red	73589-00534-00	73589-00535-00	73589-00536-00	73589-00537-00	73589-00538-00	73589-00539-00

JOBST® UltraSheer Aurora Maternity



Colour	1	2	3	4	5	6
Class 1 (18-21 mmHg)						
Cranberry	73581-00252-00	73581-00253-00	73581-00254-00	73581-00255-00	73581-00256-00	73581-00257-00
Navy Blue	73581-00258-00	73581-00259-00	73581-00260-00	73581-00261-00	73581-00262-00	73581-00263-00
Lavender	73581-00264-00	73581-00265-00	73581-00266-00	73581-00267-00	73581-00268-00	73581-00269-00
Taupe	73581-00288-00	73581-00289-00	73581-00290-00	73581-00291-00	73581-00292-00	73581-00293-00
Crystal Grey	73581-00300-00	73581-00301-00	73581-00302-00	73581-00303-00	73581-00304-00	73581-00305-00
Cherry Red	73581-00312-00	73581-00313-00	73581-00314-00	73581-00315-00	73581-00316-00	73581-00317-00
Class 2 (22-33 mmHg)						
Cranberry	73589-00450-00	73589-00451-00	73589-00452-00	73589-00453-00	73589-00454-00	73589-00455-00
Navy Blue	73589-00456-00	73589-00457-00	73589-00458-00	73589-00459-00	73589-00460-00	73589-00461-00
Lavender	73589-00462-00	73589-00463-00	73589-00464-00	73589-00465-00	73589-00466-00	73589-00467-00
Taupe	73589-00504-00	73589-00505-00	73589-00506-00	73589-00507-00	73589-00508-00	73589-00509-00
Crystal Grey	73589-00528-00	73589-00529-00	73589-00530-00	73589-00531-00	73589-00532-00	73589-00533-00
Cherry Red	73589-00546-00	73589-00547-00	73589-00548-00	73589-00549-00	73589-00550-00	73589-00551-00

JOBST® OPAQUE

HIGH CLASS COMFORT

- Smooth, lightweight design for an increased comfort
- Wicking properties keep skin breathing and dry
- Easy to don thanks to its specially integrated silicone blend softener

COMFORT AND DURABILITY

- Reinforced heel and toe for greater durability
- Superior quality yarns result in a perfectly soft and durable texture

STYLES



Knee high



Thigh high



Pantyhose



TOE VERSIONS



Open toe



Closed toe

COLOURS



Natural



Honey



Black



Espresso

BAND VARIATIONS



Silicone Dot



COMPRESSION CLASSES

15-20 mmHg

20-30 mmHg

30-40 mmHg

Not all colours and styles may be available in all compression ranges; please refer to product charts.

JOBST® Opaque Knee High



Colour	Style/Band	Toe Style	S	M	L	XL	LFC	XLFC
15-20 mmHg								
● Black		Closed	7518101	7518201	7518301	7518401		
	Petite	Closed	7518123	7518223	7518323	7518423		
	SoftFit™	Closed	7518130	7518129	7518128	7518127		
		Open	7518103	7518203	7518303	7518403		
	SoftFit™	Open	7518140	7518141	7518142	7518143		
● Natural		Closed	7518100	7518200	7518300	7518400		
	Petite	Closed	7518122	7518222	7518322	7518422		
	SoftFit™	Closed	7518126	7518125	7518124	7518131		
		Open	7518102	7518202	7518302	7518402		
	Petite	Open	7518120	7518220	7518320	7518420		
● Honey	SoftFit™	Open	7518144	7518145	7518146	7518147		
		Closed	7518109	7518209	7518309	7518409		
● Espresso	SoftFit™	Closed	7518132	7518133	7518134	7518135		
		Closed	7518110	7518210	7518310	7518410		
	SoftFit™	Closed	7518136	7518137	7518138	7518139		
20-30 mmHg								
● Black		Closed	7521703	7521803	7521903	7522003	7848923	7848823
	Petite	Closed	7521723	7521823	7521923	7522023		
	SoftFit™	Closed	7521749	7521748	7521747	7521746	7521745	7521744
		Open	7521705	7521805	7521905	7522005	7848421	7848721
	SoftFit™	Open	7521761	7521760	7521759	7521758	7521757	7521756
● Natural		Closed	7521700	7521800	7521900	7522000	7848922	7848822
	Petite	Closed	7521722	7521822	7521922	7522022		
	SoftFit™	Closed	7521743	7521742	7521741	7521740	7521739	7521738
		Open	7521704	7521804	7521904	7522004	7848420	7848720
	Petite	Open	7521720	7521820	7521920	7522020		
● Honey	SoftFit™	Open	7521755	7521754	7521753	7521752	7521750	7521751
		Closed	7521711	7521811	7521911	7522011		
● Espresso	SoftFit™	Closed	7521737	7521736	7521735	7521734		
		Closed	7521713	7521813	7521913	7522013		
	SoftFit™	Closed	7521733	7521732	7521731	7521730		
30-40 mmHg								
● Black		Closed	7522902	7523002	7523102	7523202	7848723	7848423
	Petite	Closed	7522923	7523023	7523123	7523223		
	SoftFit™	Closed	7522949	7522948	7522947	7522946	7522945	7522944
		Open	7522904	7523004	7523104	7523204	7848821	7848921
	SoftFit™	Open	7522961	7522960	7522959	7522958	7522957	7522956
● Natural		Closed	7522900	7523000	7523100	7523200	7848722	7848422
	Petite	Closed	7522922	7523022	7523122	7523222		
	SoftFit™	Closed	7522943	7522942	7522941	7522940	7522939	7522938
		Open	7522903	7523003	7523103	7523203	7848820	7848920
	Petite	Open	7522920	7523020	7523120	7523220		
● Honey	SoftFit™	Open	7522955	7522954	7522953	7522952	7522950	7522951
		Closed	7522908	7523008	7523108	7523208		
● Espresso	SoftFit™	Closed	7522937	7522936	7522935	7522934		
		Closed	7522910	7523010	7523110	7523210		
	SoftFit™	Closed	7522933	7522932	7522931	7522930		

JOBST® Opaque Thigh High



Colour	Style/Band	Toe Style	S	M	L	XL
15-20 mmHg						
● Black	Dot	Closed	7518501	7518601	7518701	7518801
	Sensitive	Closed	7768904	7768905	7768906	7768907
	Dot	Open	7518503	7518603	7518703	7518803
● Natural	Dot	Closed	7518500	7518600	7518700	7518800
	Petite/Dot	Closed	7518525	7518625	7518725	7518825
	Sensitive	Closed	7768912	7768913	7768914	7768915
	Petite/Sensitive	Closed	7768920	7768921	7768922	7768923
	Dot	Open	7518502	7518602	7518702	7518802
	Petite/Dot	Open	7518524	7518624	7518724	7518824
● Honey	Dot	Closed	7518507	7518607	7518707	7518807
	Dot	Open	7518508	7518608	7518708	7518808
● Espresso	Dot	Closed	7518509	7518609	7518709	7518809
20-30 mmHg						
● Black	Dot	Closed	7522101	7522201	7522301	7522401
	Sensitive	Closed	7769004	7769005	7769006	7769007
	Dot	Open	7522103	7522203	7522303	7522403
● Natural	Dot	Closed	7522100	7522200	7522300	7522400
	Petite/Dot	Closed	7522125	7522225	7522325	7522425
	Sensitive	Closed	7769012	7769013	7769014	7769015
	Petite/Sensitive	Closed	7769020	7769021	7769022	7769023
	Dot	Open	7522102	7522202	7522302	7522402
	Petite/Dot	Open	7522124	7522224	7522324	7522424
● Honey	Dot	Closed	7522107	7522207	7522307	7522407
	Dot	Open	7522108	7522208	7522308	7522408
● Espresso	Dot	Closed	7522109	7522209	7522309	7522409
30-40 mmHg						
● Black	Dot	Closed	7523301	7523401	7523501	7523601
	Sensitive	Closed	7769104	7769105	7769106	7769107
	Dot	Open	7523303	7523403	7523503	7523603
● Natural	Dot	Closed	7523300	7523400	7523500	7523600
	Petite/Dot	Closed	7523325	7523425	7523525	7523625
	Sensitive	Closed	7769112	7769113	7769114	7769115
	Petite/Sensitive	Closed	7769120	7769121	7769122	7769123
	Dot	Open	7523302	7523402	7523502	7523602
	Petite/Dot	Open	7523324	7523424	7523524	7523624
● Honey	Dot	Closed	7523307	7523407	7523507	7523607
	Dot	Open	7523308	7523408	7523508	7523608
● Espresso	Dot	Closed	7523309	7523409	7523509	7523609

JOBST® Opaque Pantyhose



Colour	Style/Band	Toe Style	S	M	L	XL
15-20 mmHg						
● Black		Closed	7518901	7519001	7519101	7519201
● Natural		Closed	7518900	7519000	7519100	7519200
20-30 mmHg						
● Black		Closed	7522503	7522603	7522703	7522803
● Natural		Closed	7522500	7522600	7522700	7522800
30-40 mmHg						
● Black		Closed	7523702	7523802	7523902	7524002
● Natural		Closed	7523700	7523800	7523900	7524000



JOBST® SOSOFT

CASUAL LOOK AT ITS BEST

- Stylish patterns give an edge to its laid back, contemporary design

SOFT & EASY TO DON

- Softness of yarns ensures high patient compliance
- Superior moisture management for dry and comfortable feet

STYLE



Knee high

TOE VERSION



Closed toe



COLOURS



Sand



Black



Brocade



Ribbed

VARIATIONS

JOBST® soSoft Knee High



Colour	Style	Toe Style	S	M	L	XL
15-20 mmHg						
Black	Ribbed	Closed	7456826	7456827	7456828	7456829
	Brocade	Closed	7456818	7456819	7456820	7456821
Sand	Ribbed	Closed	7456830	7456831	7456832	7456833
	Brocade	Closed	7456822	7456823	7456824	7456825
20-30 mmHg						
Black	Ribbed	Closed	7456842	7456843	7456844	7456845
	Brocade	Closed	7456834	7456835	7456836	7456837
Sand	Ribbed	Closed	7456846	7456847	7456848	7456849
	Brocade	Closed	7456838	7456839	7456840	7456841
30-40 mmHg						
Black	Ribbed	Closed	7457750	7457850	7457950	7458050
	Brocade	Closed	7456950	7457050	7457150	7457250
Sand	Ribbed	Closed	7458160	7458260	7458360	7458460
	Brocade	Closed	7457360	7457460	7457560	7457660

COMPRESSION CLASSES

15-20 mmHg

20-30 mmHg

30-40 mmHg

JOBST® SPORT

DESIGNED FOR ATHLETIC PEOPLE

- Longwearing stocking that is soft, lightweight and breathable
- Stylish and athletic design
- Odor and moisture control

STYLE



Knee high

TOE VERSION



Closed toe

COLOURS



COMING SOON

JOBST® Sport Knee High



Colour	Toe Style	S	M	L	XL
15-20 mmHg					
Dark/Light Grey	Closed	7528930	7528931	7528932	7528933
Light Grey/White	Closed	7528920	7528921	7528922	7528923
20-30 mmHg					
Dark/Light Grey	Closed	7529030	7529031	7529032	7529033
Light Grey/White	Closed	7529020	7529021	7529022	7529023

COMPRESSION CLASSES

15-20 mmHg

20-30 mmHg

JOBST® ACTIVEWEAR

IT'S ALL ABOUT COMFORT

- Extra cushioned foot with a reinforced heel and toe enhance wearing comfort and resists abrasion
- Feels like cotton with Dri-Release fabric blend
- Seamless toe

STYLE



Knee high



TOE VERSION



Closed toe



COLOURS



JOBST® ActiveWear Knee High



Colour	Toe Style	S	M	L	XL	LFC	XLFC
15-20 mmHg							
Black	Closed	110483	110484	110485	110486	110531	110532
White	Closed	110479	110480	110481	110482	110525	110526
20-30 mmHg							
Black	Closed	110493	110494	110495	110496	110533	110534
White	Closed	110489	110490	110491	110492	110527	110528
30-40 mmHg							
Black	Closed	7541103	7541203	7541303	7541403	110535	110536
White	Closed	7541102	7541202	7541302	7541402	110529	110530

COMPRESSION CLASSES

15-20 mmHg

20-30 mmHg

30-40 mmHg

JOBST® TRAVEL SOCKS

STYLISH COMPRESSION

- Ribbed pattern complements both business and casual attire

COMFORTABLE TRIP

- Soft and breathable for optimal wearing comfort
- Helps prevent leg swelling and reduces the likelihood of blood clot formation in a DVT

STYLE



Knee high

TOE VERSION



Closed toe

COLOURS



Black



Navy

JOBST® Travel Knee High



Colour	Toe Style	1	2	3	4	5
15-20 mmHg						
● Black	Closed	7884407	7884507	7884607	7884707	7927207
● Navy	Closed	7884406	7884506	7884606	7884706	7927206

COMPRESSION CLASS

15-20 mmHg

JOBST® RELIEF

EFFECTIVE & AFFORDABLE COMPRESSION

- Affordable price without compromising effectiveness of compression
- Reinforced heel and toe resist wear and tear

STYLES



Knee high



Thigh high



Pantyhose



Chap style

TOE VERSIONS



Open toe



Closed toe

COLOURS



Beige



Black

Not all colours and styles may be available in all compression ranges; please refer to product charts.

JOBST® Relief Knee High



Colour	Toe Style	S	M	L	XL	LFC	XLFC
15-20 mmHg							
● Black	Closed	114812	114813	114814	114815	114816	114817
● Beige	Closed	114806	114807	114808	114809	114810	114811
	Open	114800	114801	114802	114803	114804	114805
20-30 mmHg							
● Black	Closed	114730	114731	114732	114733	114734	114735
● Beige	Closed	114620	114621	114622	114623	114698	114624
	Open	114625	114626	114627	114628	114696	114629
30-40 mmHg							
● Black	Closed	114736	114737	114738	114739	114740	114741
● Beige	Closed	114630	114631	114632	114633	114699	114634
	Open	114635	114636	114637	114638	114697	114639

COMPRESSION CLASSES

15-20 mmHg

20-30 mmHg

30-40 mmHg

JOBST® Relief Thigh High



Colour	Band	Toe Style	S	M	L	XL
15-20 mmHg						
Black	Dot	Closed	114826	114827	114828	114829
Beige	Dot	Closed	114822	114823	114824	114825
	Dot	Open	114818	114819	114820	114821
20-30 mmHg						
Black	Dot	Closed	114212	114213	114214	114215
Beige	Dot	Closed	114208	114209	114210	114211
	No band	Closed	114640	114641	114642	114643
	Dot	Open	114200	114201	114202	114203
	No band	Open	114644	114645	114646	114647
30-40 mmHg						
Black	Dot	Closed	114220	114221	114222	114223
Beige	Dot	Closed	114216	114217	114218	114219
	No band	Closed	114648	114649	114650	114651
	Dot	Open	114204	114205	114206	114207
	No band	Open	114652	114653	114654	114655

JOBST® Relief Pantyhose



Colour	Toe Style	S	M	L	XL
20-30 mmHg					
Beige	Closed	114656	114657	114658	114659
	Open	114660	114661	114662	114663
30-40 mmHg					
Beige	Closed	114664	114665	114666	114667
	Open	114668	114669	114670	114671

JOBST® Relief Chap Style



Colour	Style	Toe Style	S	M	L	XL
20-30 mmHg						
Beige	Left Leg	Open	114672	114673	114674	114675
	Right Leg	Open	114676	114677	114678	114679
	2 Legs	Open	114680	114681	114682	114683
30-40 mmHg						
Beige	Left Leg	Open	114784	114785	114786	114787
	Right Leg	Open	114788	114789	114790	114791
	2 Legs	Open	114792	114793	114794	114795

JOBST® forMEN AMBITION

DESIGNED FOR MEN

- New heel pocket designed especially for men's anatomy for a perfect fit
- Relief zone with flat toe seam reduces pressure points for all day wearing comfort
- Soft yarns and high quality cotton for easy donning and comfortable wearing
- Odor and moisture control with C² effect

STYLE



Knee high




TOE VERSION



Closed toe

COLOURS



Black




Khaki



Brown



Navy



Grey

BAND VARIATION



Not all colours and styles may be available in all compression ranges; please refer to product charts.

COMPRESSION CLASSES

15-20 mmHg

20-30 mmHg

30-40 mmHg

JOBST® forMen Ambition Knee High



Colour	Style/Band	Toe Style	1	2	3	4	5	6
15-20 mmHg								
● Black	Regular	Closed	7766500	7766501	7766502	7766503	7766504	7766505
	Long	Closed	7766600	7766601	7766602	7766603	7766604	7766605
	Regular/SoftFit™	Closed	7764450	7764451	7764452	7764453	7764454	7764455
	Long/SoftFit™	Closed	7764750	7764751	7764752	7764753	7764754	7764755
● Khaki	Regular	Closed	7766510	7766511	7766512	7766513	7766514	7766515
	Long	Closed	7766610	7766611	7766612	7766613	7766614	7766615
	Regular/SoftFit™	Closed	7764460	7764461	7764462	7764463	7764464	7764465
	Long/SoftFit™	Closed	7764760	7764761	7764762	7764763	7764764	7764765
● Brown	Regular	Closed	7766520	7766521	7766522	7766523	7766524	7766525
	Long	Closed	7766620	7766621	7766622	7766623	7766624	7766625
	Regular/SoftFit™	Closed	7764470	7764471	7764472	7764473	7764474	7764475
	Long/SoftFit™	Closed	7764770	7764771	7764772	7764773	7764774	7764775
● Navy	Regular	Closed	7766530	7766531	7766532	7766533	7766534	7766535
	Long	Closed	7766630	7766631	7766632	7766633	7766634	7766635
	Regular/SoftFit™	Closed	7764480	7764481	7764482	7764483	7764484	7764485
	Long/SoftFit™	Closed	7764780	7764781	7764782	7764783	7764784	7764785
● Grey	Regular	Closed	7766540	7766541	7766542	7766543	7766544	7766545
	Long	Closed	7766640	7766641	7766642	7766643	7766644	7766645
	Regular/SoftFit™	Closed	7764490	7764491	7764492	7764493	7764494	7764495
	Long/SoftFit™	Closed	7764790	7764791	7764792	7764793	7764794	7764795
20-30 mmHg								
● Black	Regular	Closed	7766700	7766701	7766702	7766703	7766704	7766705
	Long	Closed	7766800	7766801	7766802	7766803	7766804	7766805
	Regular/SoftFit™	Closed	7764550	7764551	7764552	7764553	7764554	7764555
	Long/SoftFit™	Closed	7764850	7764851	7764852	7764853	7764854	7764855
● Khaki	Regular	Closed	7766710	7766711	7766712	7766713	7766714	7766715
	Long	Closed	7766810	7766811	7766812	7766813	7766814	7766815
	Regular/SoftFit™	Closed	7764560	7764561	7764562	7764563	7764564	7764565
	Long/SoftFit™	Closed	7764860	7764861	7764862	7764863	7764864	7764865
● Brown	Regular	Closed	7766720	7766721	7766722	7766723	7766724	7766725
	Long	Closed	7766820	7766821	7766822	7766823	7766824	7766825
	Regular/SoftFit™	Closed	7764570	7764571	7764572	7764573	7764574	7764575
	Long/SoftFit™	Closed	7764870	7764871	7764872	7764873	7764874	7764875
● Navy	Regular	Closed	7766730	7766731	7766732	7766733	7766734	7766735
	Long	Closed	7766830	7766831	7766832	7766833	7766834	7766835
	Regular/SoftFit™	Closed	7764580	7764581	7764582	7764583	7764584	7764585
	Long/SoftFit™	Closed	7764880	7764881	7764882	7764883	7764884	7764885
● Grey	Regular	Closed	7766740	7766741	7766742	7766743	7766744	7766745
	Long	Closed	7766840	7766841	7766842	7766843	7766844	7766845
	Regular/SoftFit™	Closed	7764590	7764591	7764592	7764593	7764594	7764595
	Long/SoftFit™	Closed	7764890	7764891	7764892	7764893	7764894	7764895

JOBST® forMen Ambition Knee High



Colour	Style/Band	Toe Style	1	2	3	4	5	6
30-40 mmHg								
● Black	Regular	Closed	7766900	7766901	7766902	7766903	7766904	7766905
	Long	Closed	7767000	7767001	7767002	7767003	7767004	7767005
	Regular/SoftFit™	Closed	7764650	7764651	7764652	7764653	7764654	7764655
	Long/SoftFit™	Closed	7764950	7764951	7764952	7764953	7764954	7764955
● Khaki	Regular	Closed	7766910	7766911	7766912	7766913	7766914	7766915
	Long	Closed	7767010	7767011	7767012	7767013	7767014	7767015
	Regular/SoftFit™	Closed	7764660	7764661	7764662	7764663	7764664	7764665
	Long/SoftFit™	Closed	7764960	7764961	7764962	7764963	7764964	7764965
● Brown	Regular	Closed	7766920	7766921	7766922	7766923	7766924	7766925
	Long	Closed	7767020	7767021	7767022	7767023	7767024	7767025
	Regular/SoftFit™	Closed	7764670	7764671	7764672	7764673	7764674	7764675
	Long/SoftFit™	Closed	7764970	7764971	7764972	7764973	7764974	7764975
● Navy	Regular	Closed	7766930	7766931	7766932	7766933	7766934	7766935
	Long	Closed	7767030	7767031	7767032	7767033	7767034	7767035
	Regular/SoftFit™	Closed	7764680	7764681	7764682	7764683	7764684	7764685
	Long/SoftFit™	Closed	7764980	7764981	7764982	7764983	7764984	7764985
● Grey	Regular	Closed	7766940	7766941	7766942	7766943	7766944	7766945
	Long	Closed	7767040	7767041	7767042	7767043	7767044	7767045
	Regular/SoftFit™	Closed	7764690	7764691	7764692	7764693	7764694	7764695
	Long/SoftFit™	Closed	7764990	7764991	7764992	7764993	7764994	7764995

JOBST® forMEN CASUAL

PERFECT LOOK FOR YOUR DAILY ACTIVITIES

- Contemporary, relaxed design that can be worn every day

ROBUST COMFORT

- Reinforced heel and toe resist wear and tear
- Superior moisture wicking properties to help keep feet dry and cool

STYLE



Knee high

TOE VERSION



Closed toe

COLOURS



Black



Khaki

JOBST® forMen Casual Knee High



Colour	Toe Style	S	M	L	XL
15-20 mmHg					
● Black	Closed	7547602	7547702	7547802	7547902
● Khaki	Closed	7547603	7547703	7547803	7547903
20-30 mmHg					
● Black	Closed	7548402	7548502	7548602	7548702
● Khaki	Closed	7548403	7548503	7548603	7548703
30-40 mmHg					
● Black	Closed	7549202	7549302	7549402	7549502
● Khaki	Closed	7549203	7549303	7549403	7549503

COMPRESSION CLASSES

15-20 mmHg

20-30 mmHg

30-40 mmHg

JOBST® forMEN

FASHIONABLE COMFORT

- Fashionable ribbed design looks like a fine men's dress sock
- Reinforced heel and toe for greater durability

STYLES



Knee high



Thigh high

TOE VERSIONS





Open toe



Closed toe

COLOURS



BAND VARIATION



Silicone Dot

Not all colours and styles may be available in all compression ranges; please refer to product charts.

COMPRESSION CLASSES

15-20 mmHg

20-30 mmHg

30-40 mmHg

JOBST® forMen Knee High



Colour	Style/Band	Toe Style	S	M	L	XL	MTall	LTall	LFC	XLFC
15-20 mmHg										
● Black		Closed	115000	115001	115002	115003		115068		
● White		Closed	115008	115009	115010	115011				
● Khaki		Closed	115012	115013	115014	115015				
20-30 mmHg										
● Black		Closed	115088	115089	115090	115091	115255	115260	115294	115295
		Open	115432	115433	115434	115435			115372	115373
● Khaki		Closed	115100	115101	115102	115103				
30-40 mmHg										
● Black		Closed	115108	115109	115110	115111	115257	115265	115296	115297
		Open	115452	115453	115454	115455			115374	115375
● Khaki		Closed	115120	115121	115122	115123				

JOBST® forMen Thigh High



Colour	Band	Toe Style	S	M	L	XL
15-20 mmHg						
● Black	Dot	Closed	115516	115517	115518	115519
● Khaki	Dot	Closed	115512	115513	115514	115515
20-30 mmHg						
● Black	Dot	Closed	115408	115409	115410	115411
● Khaki	Dot	Closed	115400	115401	115402	115403
30-40 mmHg						
● Black	Dot	Closed	115412	115413	115414	115415
● Khaki	Dot	Closed	115404	115405	115406	115407

JOBST® CUSTOM SEAMED

CUSTOMIZATION AT ITS BEST

- Available in 74 styles and options including knee, thigh, waist high, chaps, armsleeves, reinforcements and zippers to ensure patients get exactly the garment needed
- Designed for 20-65+ mmHg compression to fit a wide range of patient indications
- No latex rubber yarn to eliminate risk of allergies

STYLES

Refer to next page

COLOURS

Beige

Black

COMPRESSION CLASS

15-65+ mmHg

JOBST® Custom Seamed STYLES UPPER EXTREMITIES

Description	Code	
Glove		
Interdigital web spacer for the hand (to be worn over glove)	100536	
Mitten to wrist	100537	
Mitten to elbow	100030	
Gauntlet	100505	
Gauntlet with extended thumb	100080	
Glove to wrist	100535	
Armsleeve		
Glove to elbow	100534	
Elbow band	100016	
Half sleeve (wrist to elbow)	100515	
Half sleeve (elbow to axilla)	101155	
Half sleeve with gauntlet	100516	
Half sleeve with shoulder flap	100013	
Armsleeve (wrist to axilla)	100501	
Armsleeve with attached gauntlet	100502	
Armsleeve with attached shoulder flap	100503	
Armsleeve with gauntlet and shoulder flap	100504	
Armsleeve to axilla for stump	101140	
Armsleeve with shoulder flap for stump	101141	

Description	Code	
Mask/Chin Straps		
Face mask	100540	
Face mask, open face	101158	
Modified chin strap (extending behind the ears)	100549	
Chin strap	100550	
Vest		
Vest, sleeveless	100525	
Vest, two long sleeves	100527	
Vest, two short sleeves	100526	
Vest, one long sleeve and one short sleeve	100524	
Vest, one long sleeve, gauntlet and one short sleeve	100523	
Body Brief		
Torso band	100018	
Sleeveless body brief	100530	
Body Brief with Sleeves	100531	
Sleeveless body suit with legs	100558	
Body suit with sleeves and legs	100560	

JOBST® Custom Seamed STYLES LOWER EXTREMITIES

Description	Code	
Panty Girdle Above Knee		
Panty girdle, above knee, closed pubis	101119	
Panty girdle, above knee, open pubis	101122	
Panty Girdle Below Knee		
Panty girdle, below knee, closed pubis	101111	
Panty girdle, below knee, open pubis	101110	
Abdominal pelvic garment	101121	
Knee High		
Foot glove to ankle	100538	
Foot glove extending to knee	100539	
Ankle garment	100105	
Knee length garment	100101	
Thigh High		
Thigh band	100019	
Thigh length garment	100201	
Waist High		
Waist high, two legs, closed pubis	101101	
Waist high, two legs, open pubis	101102	
Waist high, one leg, open pubis	101103	
Waist high, one leg panty, open pubis	101112	
Waist high, one leg panty, closed pubis	101113	

Description	Code	
Waist high, one leg brief	101114	
Chest high, two legs, pregnancy	101104	
Stump garment, above knee	101130	
Stump garment, below knee	101131	
Waist high, one leg and stump garment, above or below knee, closed pubis	101134	
Waist high, one leg and stump garment, above or below knee, open pubis	101135	
Stump garment, one stump, waist high, above or below knee	101132	
Waist high, two stumps, above or below knee, closed pubis	101133	
Waist high, two stumps, above or below knee, open pubis	101136	
Waist high, one stump, above or below knee and panty	100039	
One leg, chap style garment	100035	
Two legs, chap style garment	100036	
One Stump, above or below knee, chap style garment	100037	
Two stumps, above or below knee, chap style garment	100038	

JOBST® Custom Seamed OPTIONS

Description	Code	Knee High A	Thigh High B	Waist High C	Armsleeve D	Glove E	Panty Girdle Above Knee F	Panty Girdle Below Knee G	Body Brief/ Suits H	Vest I	Mask/Chin Straps J
Gauntlet attached to body suit/brief	101170								√		
VELCRO® tabs (set of four for vest attachments)	101163			√			√	√			
Elbow lining (inner aspect only)	101167				√					√	
Elbow lining (full)	101168				√					√	
Slant inserts included in glove	101169					√					
Adjustable shoulder flap	101172				√						
Lining behind the knee (per leg)	100040		√	√				√	√		
Reinforced knee (per leg)	101186		√	√				√	√		
Extra garters (8" long - 20.3 cm)	110007										
Pull tabs (per leg)	100010	√									
Garter loops (per leg)	100044		√								
Reinforced for suspenders	101157			√			√	√			
Attached/detached suspenders	101162			√			√	√			
Reinforced top of thigh high (per leg)	101184		√								
Reinforced inner thigh and perineum	101185			√			√	√	√		
Double back panel for waist high	100022			√			√	√			
Reinforced top for waist high	100024			√			√	√			
Additional charge for circumference 50-59 $\frac{7}{8}$ " (127-151 cm)	101177			√			√	√	√	√	
Additional charge for circumference 60-69 $\frac{7}{8}$ " (152-177 cm)	100031			√			√	√	√	√	
Additional charge for circumference 70-79 $\frac{7}{8}$ " (178-201 cm)	100042			√			√	√	√	√	
2" Elastic added to custom garments	101117			√			√	√		√	
Reduced pressure abdominal panel	101161			√			√	√	√		
Double abdominal panel	101175			√			√	√	√		
Abdominal uplift panel (extending from abdomen to back)	101176			√			√	√	√		
Zipper (each opening)	101164	√	√	√	√	√	√	√			
Self-material enclosed toe	101159	√	√	√							
Soft material enclosed toe	101160	√	√	√							
Reinforced heel (per leg)	101187	√	√	√							
Reinforced heel with ankle lining (per leg)	101188	√	√	√							
Foam padding for armsleeve (3" x 5" - 2.62 x 12.7 cm)	100009				√				√	√	
Pocket for padding splint	100027	√	√	√	√	√	√	√	√	√	
Contracture seam	100176	√	√	√	√	√			√	√	
Reinforced/lined palm or dorsum	100021	√	√	√	√	√			√	√	
Lining variation	101178				√	√					
1" Dotted silicone band	101118	√	√	√	√	√			√	√	
Self-material gusset with lining	101183			√					√		
Ostomy/trachea opening	101116			√			√	√	√	√	√
Ear flap attached to mask or chin strap	100014										√
Eye flap attached to mask or chin strap	100029										√
Nose covering included in mask	101165										√
Lip covering attached to mask or chin strap	101166										√
2" Dotted silicone band	100160	√	√	√	√	√			√	√	
Fee tape transposition	100559	√	√	√	√	√	√	√	√	√	√
Seamless heel	101088	√	√	√		√					
Plastic loop	101108	√	√	√	√	√	√	√	√	√	√
Stocking slip	111345	√	√	√	√						
Zipper pull tabs	122223	√	√	√	√	√	√	√	√	√	√
Hook	-	√	√	√	√				√	√	
Cloth zipper pull	-	√	√	√	√				√	√	



JOBST® CUSTOM SEAMLESS SOFT

CIRCULAR KNIT FOR PRECISE MANAGEMENT OF VASCULAR INDICATIONS

- Controlled gradient compression for effective long-term management of mild and moderate swelling of lower extremities
- Attractive, lightweight fabric to create an anatomically fitting compression garment
- Camouflages skin imperfections

HIGH QUALITY FABRIC FOR INCREASED COMFORT

- Multifiber yarns provide unparalleled softness and high durability
- Stretchable elastic yarns for easy donning
- High air permeability allows skin to breathe

STYLES



Knee high




Thigh high



Waist high



Maternity



Chap style

COLOURS




Sand



Amber



Sun Bronze



Navy



Black

TOE VERSIONS



Open toe



Closed toe

BAND VARIATIONS



Silicone Lace



Silicone Dot



Strips

COMPRESSION CLASSES

15-20 mmHg

20-30 mmHg

30-40 mmHg

JOBST® Custom Seamless Soft

Style	Toe Style	Code
Knee high		
Knee high	Open	711020
Knee high	Closed	711040
Thigh high		
Thigh high	Open	713020
Thigh high	Closed	713040
Waist high		
Two legs	Open	715520
Two legs	Closed	715540
With fly	Open	716320
With fly	Closed	716340
Full compression	Open	718520
Full compression	Closed	718540
One leg panty	Open	719520
One leg panty	Closed	719540
Maternity	Open	717520
Maternity	Closed	717540
Chap style, one leg	Open	714020
Chap style, one leg	Closed	714040
Chap style, two legs	Open	715320
Chap style, two legs	Closed	715340

OPTIONS

Description	Code
Sand	899115
Amber	899119
Sun Bronze	899111
Navy	899113
Black	899114
5 cm dotted silicone band (Thigh/Knee)	799102
6 cm lace silicone band (Thigh)	799103
6 cm strips silicone band (Thigh)	799104
2.5 cm dotted silicone band (Knee)	799105
Open pubis (Waist)	799107
5 cm waistband (Panty)	799109
2.5 cm waistband (Panty)	799115

JOBST® BELLAVER®

AFFORDABLE MANAGEMENT OF MODERATE TO SEVERE LOWER EXTREMITY SWELLING

- Custom garment with firm fabric to provide controlled stretch and resistance
- Easy measuring system requires up to 30% fewer measurements than a conventional seamed garment

STYLES



Knee high




Thigh high



Waist high



Maternity



Chap style

COLOURS



Sand



Sun Bronze



Black

TOE VERSIONS



Open toe



Closed toe

BAND VARIATIONS



Silicone Lace



Silicone Dot



Strips

COMPRESSION CLASSES

20-30 mmHg

30-40 mmHg

JOBST® Bellavar®

Style	Toe Style	Code
Knee high		
Knee high	Open	711021
Knee high	Closed	711041
Thigh high		
Thigh high	Open	713021
Thigh high	Closed	713041
Waist high		
Two legs	Open	715521
Two legs	Closed	715541
With fly	Open	716322
With fly	Closed	716342
Full compression	Open	718523
Full compression	Closed	718543
One leg panty	Open	719521
One leg panty	Closed	719541
One leg panty, with fly	Open	716321
One leg panty, with fly	Closed	716341
Maternity	Open	717521
Maternity	Closed	717541
Chap style, one leg	Open	714021
Chap style, one leg	Closed	714041

OPTIONS

Description	Code
Sand	899115
Sun Bronze	899111
Black	899114
5 cm dotted silicone band (Thigh/Knee)	799102
6 cm lace silicone band (Thigh)	799103
6 cm strips silicone band (Thigh)	799104
2.5 cm dotted silicone band (Knee)	799105
Open pubis (Waist)	799107

JOBST® ELVAREX®

EFFECTIVE DAYTIME THERAPY

- Custom manufacturing to ensure a contoured fit to the shape of the limb
- Controlled gradient compression for the effective long-term management of lymphedema

FLAT KNIT COMFORT

- Flat knitted construction that allows greater air permeability
- Excellent durability for extended wear

STYLES




COLOURS

Upper Extremity



Lower Extremity



COMPRESSION CLASSES

Upper Extremity Class 1, 2 and 2F

Lower Extremity Class 1, 2, 3, 3F, 4 and 4S

OPTIONS Refer to page 44

Not all colours and styles may be available in all compression ranges; please refer to product charts.



JOBST® ELVAREX® SOFT

COMFORT & SOFTNESS PATIENTS DESERVE

- For the management of mild to moderate edema
- New knitting construction and softer yarn provides easy donning
- Soft surface that is extremely comfortable to wear
- Comfortable edging seam provides the required stability and durability of the fabric edge

STYLES

Armsleeve

Gauntlet

Glove

Glove with extension piece

Knee high with slant open toe

Thigh high with slant open toe

2 legs pantyhose

Foot cap

Capri

Biker short

Chap style slant open toe

COLOURS

Upper Extremity

Black

Beige

Cranberry

Honey

Lower Extremity

Black

Beige

COMPRESSION CLASSES

- Upper Extremity Class 1 and 2
- Lower Extremity Class 1, 2 and 3

OPTIONS Refer to page 44

Not all colours and styles may be available in all compression ranges; please refer to product charts.

JOBST® ELVAREX® PLUS

ADVANCED TECHNOLOGY FOR INCREASED COMFORT

- 3D knitting technique reproducing each body shape accurately for more freedom of movement
- Seamless material for consistent flexibility and maximum range of motion
- Superior wearing comfort without pressure points
- Same JOBST® Elvarex® knitted fabric providing excellent breathability and micro-massage effect
- Ideal for patients with sensitive skin or living an active lifestyle

STYLES



Gauntlet



Glove



Foot cap

COLOURS


Foot



Black




Beige



Navy



Cocoa



Grey

Hand



Black



Beige



Caramel



Cranberry



Honey

Not all colours and styles may be available in all compression ranges; please refer to product charts.

COMPRESSION CLASSES

Upper Extremity **Class 1 and 2**

JOBST® RELAX

REVOLUTIONARY NIGHTTIME COMPRESSION THERAPY

- Customized, flat-knitting technology provides consistent gradient compression and stays in place
- The textured surface is designed to create differentiated areas of pressure that stimulate lymph flow
- The integrated zipper option with stay-down puller facilitates donning and doffing

FOR A GOOD NIGHT'S SLEEP

- The monofilament spacer layer provides cushioning to guarantee comfort during the night
- Coolmax® yarns and advanced knitting techniques increase breathability to keep patients cooler and more comfortable

STYLES



Armsleeve



Armsleeve
w/ gauntlet
attached



Knee high with
slant open toe



Thigh high with
slant open toe

COLOURS

Beige

Rosé

OPTIONS Refer to page 44

COMPRESSION CLASSES

Upper Extremity 15-20 mmHg

Lower Extremity 15-20 and 20-30 mmHg

LYMPHOLOGY CUSTOM GARMENTS OPTIONS

Description		JOBST® Elvarex®	JOBST® Elvarex® Soft	JOBST® Relax
Options				
Shoulder cap, strap		√		
Shoulder cap, bra loop		√		
Comfort zone		√		
Pocket		√	√	
Silicone band		√	√	
SoftFit™ band 		√		
Fly for men		√		
Open pubis		√		
Ajustable waist band		√	√	
Slanted closed toe		√	√	
Straight closed toe		√	√	
Open toe		√	√	
Zipper		√		√
Profile		√		
T-Heel		√	√	



JOBST® BELLA™ LITE

EFFECTIVE PREVENTION/MANAGEMENT OF MILD EDEMA

- Soft and air permeable
- Lightweight, knit construction
- Reduced compression at the wrist provides an effective and comfortable fit when the armsleeve and gauntlet are worn together
- Easy to put on, comfortable to wear

STYLES



Armsleeve



Gauntlet




Silicone Lace

BAND VARIATION

COLOURS



Beige



Black

JOBST® Bella™ Lite Armsleeve



Colour	Style	Band	S	M	L
15-20 mmHg					
	Regular		7527302	7527402	7527502
	Long		101328	101329	101330
	Regular	Silicone	7527303	7527403	7527503
	Long	Silicone	101331	101332	101333
20-30 mmHg					
	Regular		7527605	7527705	7527805
	Long		101313	101314	101315
	Regular	Silicone	7527606	7527706	7527806
	Long	Silicone	101334	101335	101336
	Regular		101340	101341	101342

COMPRESSION CLASSES

15-20 mmHg

20-30 mmHg

SIZING CHART

Measurements (cm)	S	M	L
Base of fingers	15-19	19-22	21-26
Wrist	14-20	16-22	18-24
Widest forearm	18-23	23-29	29-34
Widest bicep	22-40	27-45	33-51
Regular	≤ 43		
Long	> 43		

JOBST® Bella™ Lite Gauntlet



Colour	S	M	L
15-20 mmHg			
	101325	101326	101327
20-30 mmHg			
	101319	101320	101321
	101358	101359	101360

JOBST® BELLA™ STRONG

THERAPEUTICALLY EFFECTIVE FOR MODERATE EDEMA

- Covers unrivalled range of patient dimensions with 10 sizes per compression class
- Easy to don, armsleeve comes with donning device
- Flat wrist seam reduces rolling
- Cotton-like feel for more comfort and better compliance

STYLES



Armsleeve



Gauntlet



COLOURS

Natural

BAND VARIATION



Silicone Lace

COMPRESSION CLASSES

15-20 mmHg

20-30 mmHg

30-40 mmHg

JOBST® Bella™ Strong Armsleeve



Colour	Band	1	2	3	4	5	6	7	8	9	10
15-20 mmHg											
Natural		7272711	7272712	7272713	7272714	7272715	7272716	7272717	7272718	7272719	7272720
	Silicone	7272701	7272702	7272703	7272704	7272705	7272706	7272707	7272708	7272709	7272710
20-30 mmHg											
Natural		7273611	7273612	7273613	7273614	7273615	7273616	7273617	7273618	7273619	7273620
	Silicone	7273601	7273602	7273603	7273604	7273605	7273606	7273607	7273608	7273609	7273610
30-40 mmHg											
Natural		7277411	7277412	7277413	7277414	7277415	7277416	7277417	7277418	7277419	7277420
	Silicone	7277401	7277402	7277403	7277404	7277405	7277406	7277407	7277408	7277409	7277410

JOBST® Bella™ Strong Gauntlet



Colour	1	2	3	4	5	6	7	8
15-20 mmHg								
Natural	7277700	7277701	7277702	7277703	7277704	7277705	7277706	7277707
20-30 mmHg								
Natural	7369440	7369441	7369442	7369443	7369444	7369445	7369446	7369447

SIZING CHART

Measurements (cm)	1	2	3	4	5	6	7	8	9	10
Gauntlet										
Base of fingers	15-17	17-19	17-19	19-21	19-21	21-23	21-23	23-26		
Wrist	14-18	14-18	18-21	14-18	18-21	16-19	19-22	19-22		≤ 43
Armsleeve										
Wrist	14-18	14-18	14-18	16-19	17-20	14-18	17-20	18-21	16-19	19-22
Widest forearm	17-21	17-21	21-25	21-25	21-25	25-29	25-29	25-29	29-33	29-33
Widest bicep	22-31	31-41	25-34	28-38	34-43	29-38	33-43	38-48	32-42	36-46
Length	≤ 43									

JOBST® FARROWWRAP®

SHORT STRETCH COMPRESSION WRAP

- Helps prevent limb from swelling during the day
- Easy to don for patients with limited mobility
- Easy to use adjustable bands
- Supports patients' independence

STYLES



Footpiece



Legpiece



Thighpiece



Armpiece



Gauntlet

COLOURS



Black



Tan

VARIATIONS

- Ready-to-Wear (RTW)
- Custom-Made (CM)
- Trim-to-Fit (TTF)

COMPRESSION CLASSES

20-30 mmHg

30-40 mmHg

JOBST® FarrowWrap®

STYLES UPPER EXTREMITIES

Description	Colour	XS	S	M	L	XL
Armpiece						
LITE						
Short Right	Tan		7665012	7665014	7665016	
Short Left	Tan		7665011	7665013	7665015	
Regular Right	Tan		7665006	7665008	7665010	
Regular Left	Tan		7665005	7665007	7665009	
Long Right	Tan		7665000	7665002	7665004	
Long Left	Tan		7665035	7665001	7665003	
Short Right	Black		7665030	7665032	7665034	
Short Left	Black		7665029	7665031	7665033	
Regular Right	Black		7665024	7665026	7665028	
Regular Left	Black		7665023	7665025	7665027	
Long Right	Black		7665018	7665020	7665022	
Long Left	Black		7665017	7665019	7665021	
TTF Right (includes Velcro®)	Tan			7612002		7612003
TTF Left (includes Velcro®)	Tan			7612000		7612001
Extra Velcro® for TTF Armpiece	7612004					
Gauntlet						
Ambidextrous Hand Gauntlet with Foam Squares	Tan	7665300	7665301	7665302	7665303	



LITE RTW



LITE Hand Gauntlet



CLASSIC Custom

STYLES LOWER EXTREMITIES

Description	Colour	XS	S	M	L	XL
Thighpiece						
LITE						
Short	Tan	7665505	7665506	7665507	7665508	7665509
Regular	Tan	7665500	7665501	7665502	7665503	7665504
Long	Tan	7665510	7665511	7665512	7665513	7665514
CLASSIC						
Short	Tan	7665806	7665807	7665808	7665809	7665810
Regular	Tan	7665800	7665801	7665802	7665803	7665804
Long	Tan	7665812	7665813	7665814	7665815	7665816
STRONG						
Short	Tan	7664905	7664906	7664907	7664908	7664909
Regular	Tan	7664900	7664901	7664902	7664903	7664904
Long	Tan	7664910	7664911	7664912	7664913	7664914



LITE STRONG

JOBST® FarrowWrap®

STYLES LOWER EXTREMITIES

Description	Colour	XS	S	M	L	XL
Legpiece						
LITE						
Regular	Tan	7665400	7665401	7665402	7665403	7665404
Long	Tan	7665405	7665406	7665407	7665408	7665409
TTF Legpiece (includes Velcro®)	Tan			7612300		
Extra Velcro® for TTF Legpiece	Tan			7612301		
Extra Band for TTF Legpiece	Tan			7612302		
CLASSIC						
Regular	Tan	7665700	7665701	7665702	7665703	7665704
Long	Tan	7665705	7665706	7665707	7665708	7665709
STRONG						
Regular	Tan	7664800	7664801	7664802	7664803	7664804
Long	Tan	7664805	7664806	7664807	7664808	7664809
TTF Legpiece (includes Velcro®)	Tan			7611800		7611801
Extra Velcro® for TTF Legpiece	Tan			7611802		
Extra Band for TTF Legpiece	Tan			7611803		7611804
BASIC*						
Regular (includes Velcro®)	Tan	7666000	7666001	7666002	7666003	7666004
Long (includes Velcro®)	Tan	7666005	7666006	7666007	7666008	7666009
Extra Velcro®				7666010		
4000*						
Regular	Tan	7666100	7666101	7666102	7666103	
Long	Tan	7666104	7666105	7666106	7666107	
Regular	Black	7666108	7666109	7666110	7666111	
Long	Black	7666112	7666113	7666114	7666115	
Footpiece						
LITE						
Regular	Tan	7665205	7665206	7665207	7665208	7665209
Long	Tan	7665200	7665201	7665202	7665203	7665204
TTF Footpiece (includes Velcro®)	Tan			7612100		
Extra Velcro® for TTF Footpiece	Tan			7667500		
CLASSIC						
Regular	Tan	7665605	7665606	7665607	7665608	7665609
Long	Tan	7665600	7665601	7665602	7665603	7665604
STRONG						
Regular	Tan	7664604	7664605	7664606	7664607	7664608
Long	Tan	7664609	7664600	7664601	7664602	7664603
TTF Footpiece (includes Velcro®)	Tan			7611600		
Extra Velcro® for TTF Footpiece	Tan			7667500		
BASIC*						
Regular	Tan	7665905	7665906	7665907	7665908	7665909
Long	Tan	7665900	7665901	7665902	7665903	7665904



CLASSIC



4000



CLASSIC

*Free pair of JOBST® FarrowHybrid® ADI Wide Foot Compression Liners included

JOBST® FarrowWrap®

UPPER EXTREMITIES SIZING CHART

Measurements (cm)	XS	S	M	L	XL
Armpiece					
G - Axilla circumference		22-31	29-39	32-45	
E - Elbow crease circumference		20-27	25-34	30-40	
C - Wrist crease circumference		14-18	16-21	19-25	
C-G Short*	40-43				
C-G Regular*	44-47				
C-G Long*	48-51				
Ambidextrous Hand Gauntlet					
Base of fingers circumference	16-18	18-20	20-22	22-24	
Fold of thumb circumference	17-19	19-21	21-23	23-26	
Wrist circumference	15-18	17-20	19-22	22-24	
Trim-To-Fit Armpiece					
G - Axilla circumference			25-44		30-48
F - Widest bicep circumference			25-42		25-45
E - Elbow crease circumference			20-38		25-42
D - Widest forearm circumference			20-36		20-40
C - Wrist crease circumference			15-31		15-35
E-G Elbow crease to axilla length			15-23		20-28
C-E Wrist to elbow crease length			20-25		20-28

*C-G: Length from the wrist to the axilla

ACCESSORIES

Description	Colour	S	M	L	XL
JOBST® FarrowHybrid® Foot Compression					
ADI Regular 20-30 mmHg	Tan	7666400	7666401	7666402	7666403
ADI Wide 20-30 mmHg	Tan	7666412	7666413	7666414	7666415
ADI Wide 20-30 mmHg	Black	7666416	7666417	7666418	7666419
ADII Regular 30-40 mmHg	Tan	7666404	7666405	7666406	7666407
ADII Wide 30-40 mmHg	Tan	7666408	7666409	7666410	7666411
TG® Soft Liners					
Knee high open toe	Tan	7666900	7666902	7666904	
Thigh high open toe	Tan	7667000	7667002	7667004	
Arm	Tan	7667006	7667008	7667010	
Flesh Tone Silver Liners					
Knee high	Tan	7666800		7666802	
Arm	Tan	7666600		7666700	
JOBST® Farrow GarmentGrip™ Neoprene Roll					
14 cm x 80 cm	Tan		7667200		
14 cm x 120 cm	Tan		7667201		
14 cm x 550 cm	Tan		7667202		
20 cm x 120 cm	Tan		7667203		
30 cm x 150 cm	Tan		7667204		

If measurements don't fit in sizing charts, please refer to the JOBST® Custom-Made order forms.

LOWER EXTREMITIES SIZING CHART

Measurements (cm)	XS	S	M	L	XL
Thighpiece					
Groin circumference	60-70	67-75	70-80	75-90	85-100
Above knee circumference	47-52	53-59	60-66	67-73	74-80
E-G Short			16-20		
E-G Regular			21-26		
E-G Long			27-31		
Legpiece (except JOBST® FarrowWrap® 4000)					
Mid-calf circumference	36-43	42-50	48-58	53-63	58-68
Ankle circumference	21-25	25-30	30-36	36-42	42-50
A-D Regular*	33-37	35-39	37-41	39-43	39-43
A-D Long*	38-41	40-43	42-45	44-47	44-47
B-D Regular**	28	30	31	33	33
B-D Long**	33	35	36	38	38
Legpiece JOBST® FarrowWrap® 4000					
Mid-calf circumference	25-38	32-48	37-56	43-67	
Ankle circumference	17-23	19-28	25-37	29-43	
B-D Regular**			32		
B-D Long**			36		
Trim-To-Fit Legpiece					
L4 and L5 circumferences			20-60		45-90
L3 circumference			25-42		45-90
L2 circumference			20-38		35-80
B Ankle circumference			20-36		35-70
B-D Posterior leg length			15-31		27-35
Footpiece					
Midfoot circumference	22-24	25-27	28-30	31-34	35-40
Regular	16-17	17-18	19-20	20-21	22-23
Long	18-19	19-20	21-22	22-23	24-25
Trim-To-Fit Footpiece					
Forefoot circumference			20-30		
Heel to forefoot length			14-21		



*From floor to 2 fingers below knee crease

**B-D: Actual height of the garment


ACCESSORIES SIZING CHART

Measurements (cm)	S	M	L	XL
JOBST® FarrowHybrid® Foot Compression				
Men shoe size	3-6	7-9	10-12	13-15
Women shoe size	4-7	8-10	11-13	14-16
Regular mid-foot circumference			20-27	
Wide mid-foot circumference			28-35	
TG® Soft Liners				
Arm, calf or groin circumference	< 40	41-70	71-120	
Flesh Tone Silver Liners				
Knee high	19-48		48-70	
Arm	19-48		48-70	

USE THIS CHART AS A GUIDE TO HELP YOU



SELECT THE RIGHT JOBST® FARROWWRAP®



COMPRILAN®

SHORT STRETCH COMPRESSION BANDAGE



Unlike long stretch bandages, Comprilan® short stretch provides high working pressure and low resting pressure, making it better tolerated by patients. It enables for cost-effective treatment with easy reapplication.

INDICATIONS

- Venous and mixed leg ulcers
- Edema
- Lymphedema
- Thrombophlebitis
- Post sclerotherapy

FEATURES AND BENEFITS

- Better tolerated by patients
- Safe for mixed ulcers ($0.6 < \text{ABPI}$)
- Cost-effective when frequent dressing changes required (washable and reusable up to 10 times)
- Easy to apply and safe to use, even with minimal training
- As effective as multi-layer bandaging

Dimensions	Code	Box Contents	Selling Unit
4 cm x 5 m	7718700	1 roll	Box
6 cm x 5 m	0102600	1 roll	Box
8 cm x 5 m	0102700	1 roll	Box
10 cm x 5 m	0102800	1 roll	Box
12 cm x 5 m	0102900	1 roll	Box
10 cm x 10 m	7718800	1 roll	Box
12 cm x 10 m	7718900	1 roll	Box

COMPRIFOAM®

FOAM BANDAGE



Comprifoam® open cell foam bandage is designed to evenly distribute the pressure from compression bandages over extremities, thus preventing constrictions and ensuring an effective compression therapy.

INDICATIONS

- Padding under compression bandages

FEATURES AND BENEFITS

- Breathable
- Non-sliding
- Tear resistant
- Easy for patients to self bandage
- Anatomically conforming allowing freedom of movement
- Resistant to body oils, lotions and perspiration
- Reusable
- May be steam sterilized

Dimensions	Code	Box Contents	Selling Unit
10 cm x 2.5 m x 0.4 cm	7529400	1 roll	Box
12 cm x 2.5 m x 0.4 cm	7529500	1 roll	Box

ARTIFLEX®

NON-WOVEN PADDING BANDAGE



Artiflex® is a multipurpose padding bandage that ensures even distribution of pressure in phlebological or lymphedema applications. Also useful as under padding of plaster or synthetic casts.

INDICATIONS

- Padding under compression bandages for sensitive skin
- Under padding of plaster and synthetic casts

FEATURES AND BENEFITS

- Non-absorbent
- Applies easily without creasing
- Provides added comfort in wrapping for lymphedema
- Breathable
- Autoclavable

Dimensions	Code	Box Contents	Selling Unit
10 cm x 3 m	09046	30 rolls	Box
15 cm x 3 m	09047	20 rolls	Box

FOAM RUBBER PAD



Dimensions	Code	Box Contents	Selling Unit
9 cm x 1 cm	78496	Unit	Box

JOBST® COMPRI2

2-LAYER BANDAGE SYSTEM



JOBST® Compri2, the new 2-layer short stretch compression system, addresses the two key elements required to ensure successful compression therapy: correct compression level easily applied and continuous patient adherence.

INDICATIONS

- Venous and mixed leg ulcers
- Lymphedema

FEATURES AND BENEFITS

- Provides sustained graduated compression of 40 mmHg (JOBST® Compri2) or 20-30 mmHg (JOBST® Compri2 Lite) for up to 7 days
- Easy to apply with 2-way stretch padding and compression indicators
- Corrective action / re-application possible
- Comfortable skin-friendly padding
- Thin system
- Online support (application videos, micro-site for patients and professionals)

Description	Code	Ankle size	Box Contents	Selling Unit
JOBST® Compri2	7627100	18-25 cm	1 kit	Box
	7627101	25-32 cm	1 kit	Box
JOBST® Compri2 lite	7627102	18-25 cm	1 kit	Box
	7627103	25-32 cm	1 kit	Box

JOBST® ULCERCARE™

READY-TO-WEAR COMPRESSION SYSTEM



JOBST® UlcerCARE™ enables effective management of a venous leg ulcer while offering ease of use at the same time – the perfect combination for providing patients with a higher quality of life.

INDICATIONS

- Venous leg ulcers

FEATURES AND BENEFITS

- Provides graduated compression of 40 mmHg
- High stiffness factor for effective compression
- Easy handling by patients or caregivers
- Contains silk in the liner for ease of donning
- Supports patients' independence
- Compatible with common footwear and clothing to enable a normal lifestyle

Style	Colour	S	M	L	XL	2XL	3XL	4XL	Box Contents	Selling Unit
No Zipper	Beige	7363021	7363022	7363023	7363024	7363025	7363026	7363027	1 stocking & 2 liners	Box
No Zipper	Black	7363041	7363042	7363043	7363044	7363045	7363046	7363047	1 stocking & 2 liners	Box
Left Zipper	Beige	7363121	7363122	7363123	7363124	7363125	7363126	7363127	1 stocking & 2 liners	Box
Left Zipper	Black	7363141	7363142	7363143	7363144	7363145	7363146	7363147	1 stocking & 2 liners	Box
Right Zipper	Beige	7363161	7363162	7363163	7363164	7363165	7363166	7363167	1 stocking & 2 liners	Box
Right Zipper	Black	7363171	7363172	7363173	7363174	7363175	7363176	7363177	1 stocking & 2 liners	Box
Replacement Liners										
	White	7363221	7363222	7363223	7363224	7363225	7363226	7363227	3 liners	Box

SIZING CHART

Measurements (cm)	S	M	L	XL	2XL	3XL	4XL
Ankle	18-21	21-24	24-27	27-30	30-33	33-36	36-39
Calf	29-35	34-41	40-47	46-54	48-56	50-58	52-60



JOBST® SURGICAL VEST

IDEAL SOLUTION FOR POST-OPERATIVE RECOVERY

- Supports breast and chest tissue after breast surgery
- Features break-resistant nylon zipper and Velcro® hook and eye closures in the front for easy adjustments
- Luxuriously soft fabric
- Minimizes hematoma formation
- Holds surgical implants in place
- Allows easy examination of the surgical site
- Helps to reduce tension on surgical incision
- Available with or without cups

SIZING CHART

Measurements (cm)	1	2	3	4	5
Chest (under arm)	79-89	89-99	99-109	109-119	120-130
Bra size	30C, 32B, 32C, 34A, 34B, 36A	34C, 34D, 36B, 36C, 38A, 38B	34DD, 36D, 36DD, 38C, 38D, 40C	38DD, 40D, 40DD, 42C, 42D, 44C	42DD, 44D, 44DD, 46C, 46D, 48C

JOBST® Compression Vest

Description	1	2	3	4	5
With cups	111901	111902	111903	111904	111905
Without cups	111906	111907	111908	111909	111910
With left cup only	111911	111912	111913	111914	111915
With right cup only	111916	111917	111918	111919	111920

IT STAYS®!

ROLL-ON BODY ADHESIVE



It Stays!® Roll-on Body Adhesive is specially designed to help keep JOBST® stockings in place. Ideal for thigh length stockings when garter belts are undesirable or silicone band is not available.

Description	Code	Box Contents	Case Contents
2 oz. (59 mL) bottle	112013	1 each	N/A
2 oz. (59 mL) bottle	112014	N/A	12 bottles

JOBST® STOCKING DONNER



This product is ideal for patients with arthritis and limited bending motion. Use with closed as well as open toe compression stockings.

Description	Code	Box Contents
JOBST® Stocking Donner	110913	1 each

JOLASTIC®

SPECIAL WASHING SOLUTION



Washing solution for elastic garments formulated to remove soil, body acids, and skin salts quickly and easily without damage to fabric. This specially designed solution will help extend the life of all elastic garments.

Description	Code	Box Contents	Case Contents
12 oz. (355 mL) bottle	112001	1 each	N/A
1 qt. (946 mL) bottle	112003	1 each	N/A
4 oz. (118 mL) bottle	130999	N/A	12 bottles

STOCKING SLIP



Reduces friction between the compression stocking and the skin allowing the compression stocking to slip on easily. Use with open toe compression stockings only.

Description	Code	Bag Contents
JOBST® Stocking Slip	112018	6 each

FOOT REST



Provides convenience when measuring for compression garments by elevating the patient's foot at the heel.

Description	Code	Box Contents	Case Contents
Foot Rest	110912	1 each	N/A

JOBST® DONNING GLOVES



JOBST® gloves create friction and provide grip to help donning and positioning evenly the fabric on the leg or the arm.

Description	Code	Bag Contents
Small	131202	1 pair
Medium	131203	1 pair

FIBER CONTENTS

	Nylon	Spandex	Silicone	Polyester	Cotton	Wool	Elastodiene
JOBST® Ultrasheer							
Knee High 15-20 mmHg	75%	25%					
Knee High 20-30 mmHg	74%	26%					
Knee High 30-40 mmHg	60%	40%					
Knee High SoftFit™ Band 15-20 mmHg	67%	27%	6%				
Knee High SoftFit™ Band 20-30 mmHg	67%	27%	6%				
Knee High SoftFit™ Band 30-40 mmHg	54%	42%	4%				
Thigh High Silicone Lace Band 15-20 mmHg	60%	25%	15%				
Thigh High Silicone Lace Band 20-30 mmHg	59%	29%	12%				
Thigh High Silicone Lace Band 30-40 mmHg	55%	37%	8%				
Thigh High Silicone Dot Band 15-20 mmHg	66%	19%	15%				
Thigh High Silicone Dot Band 20-30 mmHg	64%	24%	12%				
Thigh High Silicone Dot Band 30-40 mmHg	56%	34%	10%				
Thigh High Sensitive Band 15-20 mmHg	75%	21%	4%				
Thigh High Sensitive Band 20-30 mmHg	71%	26%	3%				
Thigh High Sensitive Band 30-40 mmHg	64%	33%	3%				
Pantyhose 15-20 mmHg	72%	28%					
Pantyhose 20-30 mmHg	73%	27%					
Pantyhose 30-40 mmHg	61%	39%					
Maternity 15-20 mmHg	72%	28%					
Maternity 20-30 mmHg	76%	24%					
Aurora Edition 15-20 mmHg	58%	42%					
Aurora Edition 20-30 mmHg	53%	47%					
JOBST® Opaque							
Knee High 15-20 mmHg	78%	22%					
Knee High 20-30 mmHg	78%	22%					
Knee High 30-40 mmHg	70%	30%					
Knee High SoftFit™ Band 15-20 mmHg	75%	21%	4%				
Knee High SoftFit™ Band 20-30 mmHg	72%	24%	4%				
Knee High SoftFit™ Band 30-40 mmHg	62%	34%	4%				
Thigh High 15-20 mmHg	72%	18%	10%				
Thigh High 20-30 mmHg	70%	16%	14%				
Thigh High 30-40 mmHg	62%	28%	10%				
Thigh High Sensitive Band 15-20 mmHg	76%	21%	3%				
Thigh High Sensitive Band 20-30 mmHg	76%	21%	3%				
Thigh High Sensitive Band 30-40 mmHg	66%	31%	3%				
Pantyhose 15-20 mmHg	76%	24%					
Pantyhose 20-30 mmHg	76%	24%					
Pantyhose 30-40 mmHg	66%	34%					
JOBST® soSoft							
15-20 mmHg	81%	19%					
20-30 mmHg	74%	26%					
30-40 mmHg	74%	26%					
JOBST® Sport							
15-20 mmHg	60%	10%		30%			
20-30 mmHg	60%	10%		30%			
JOBST® ActiveWear							
15-20 mmHg	20%	9%		60%	11%		
20-30 mmHg	20%	9%		60%	11%		
30-40 mmHg	20%	12%		58%	10%		
JOBST® Travel Socks	55%	12%			33%		

	Nylon	Spandex	Silicone	Polyester	Cotton	Wool	Elastodiene
JOBST® Relief							
Knee High Closed Toe 15-20 mmHg	80%	20%					
Knee High Closed Toe 20-30 mmHg	77%	23%					
Knee High Closed Toe 30-40 mmHg	72%	28%					
Knee High Open Toe 15-20 mmHg	80%	20%					
Knee High Open Toe 20-30 mmHg	75%	25%					
Knee High Open Toe 30-40 mmHg	70%	30%					
Thigh High Silicone Band 15-20 mmHg	74%	18%	8%				
Thigh High Silicone Band 20-30 mmHg	70%	23%	7%				
Thigh High Silicone Band 30-40 mmHg	66%	28%	6%				
Thigh High Without Band Closed Toe 20-30 mmHg	74%	26%					
Thigh High Without Band Closed Toe 30-40 mmHg	67%	33%					
Thigh High Without Band Open Toe 20-30 mmHg	73%	27%					
Thigh High Without Band Open Toe 30-40 mmHg	68%	32%					
Pantyhose Closed Toe 20-30 mmHg	75%	25%					
Pantyhose Closed Toe 30-40 mmHg	75%	25%					
Pantyhose Open Toe 20-30 mmHg	78%	22%					
Pantyhose Open Toe 30-40 mmHg	78%	22%					
Chap 20-30 mmHg	78%	22%					
Chap 30-40 mmHg	78%	22%					
JOBST® forMen Ambition							
Knee High 15-20 mmHg	63%	21%		1%	15%		
Knee High 20-30 mmHg	59%	25%		1%	15%		
Knee High 30-40 mmHg	57%	27%		1%	15%		
Knee High SoftFit™ Band 15-20 mmHg	58%	22%	2%		18%		
Knee High SoftFit™ Band 20-30 mmHg	52%	28%	2%		18%		
Knee High SoftFit™ Band 30-40 mmHg	50%	30%	2%		18%		
JOBST® forMen Casual							
15-20 mmHg	46%	14%		35%		5%	
20-30 mmHg	47%	16%		33%		4%	
30-40 mmHg	48%	16%		31%		5%	
JOBST® forMen							
Knee High Closed Toe 15-20 mmHg	76%	24%					
Knee High Closed Toe 20-30 mmHg	73%	27%					
Knee High Closed Toe 30-40 mmHg	73%	27%					
Knee High Open Toe 20-30 mmHg	71%	29%					
Knee High Open Toe 30-40 mmHg	71%	29%					
Thigh High 15-20 mmHg	73%	20%	7%				
Thigh High 20-30 mmHg	71%	25%	4%				
Thigh High 30-40 mmHg	71%	25%	4%				
JOBST® Bella™ Strong							
Armsleeve 15-20 mmHg	80%	20%					
Armsleeve 20-30 mmHg	76%	24%					
Armsleeve 30-40 mmHg	74%	26%					
Armsleeve with Silicone Band 15-20 mmHg	74%	18%	8%				
Armsleeve with Silicone Band 20-30 mmHg	72%	20%	8%				
Armsleeve with Silicone Band 30-40 mmHg	69%	23%	8%				
Gauntlet 15-20 mmHg	91%	9%					
Gauntlet 20-30 mmHg	91%	9%					
JOBST® Bella™ Lite							
Armsleeve 15-20 mmHg	84%	16%					
Armsleeve 20-30 mmHg	80%	20%					
Armsleeve with Silicone Band 15-20 mmHg	80%	17%	3%				
Armsleeve with Silicone Band 20-30 mmHg	75%	21%	3%				
Gauntlet 15-20 mmHg	84%	16%					
Gauntlet 20-30 mmHg	86%	14%					
JOBST® Custom Seamed							
LYCRA and Nylon	X	X					
JOBST® Custom Seamless Soft	X	X					
JOBST® Bellavar®	X	X					
JOBST® Elvarex®							
Upper Extremities	X	X					
Lower Extremities	X	X			X		
JOBST® Elvarex® Soft	X	X					
JOBST® Elvarex® Plus	X	X					X
JOBST® Relax	X	X		X			

BSN MEDICAL BRANDS

Artiflex®	55
Comprifoam®	55
Comprilan®	54
JOBST® ActiveWear	22
JOBST® Bella™ Lite	45
JOBST® Bella™ Strong	46
JOBST® Bellavar®	38
JOBST® Compri2	56
JOBST® Custom Seamed	32
JOBST® Custom Seamless Soft	36
JOBST® Elvarex®	40
JOBST® Elvarex® Plus	42
JOBST® Elvarex® Soft	41
JOBST® FarrowWrap®	48
JOBST® forMen	30
JOBST® forMen Ambition	26
JOBST® forMen Casual	29
JOBST® Opaque	16
JOBST® Relax	43
JOBST® Relief	24
JOBST® soSoft	20
JOBST® Sport	21
JOBST® Surgical vest	58
JOBST® Travel Socks	23
JOBST® UlcerCARE™	57
JOBST® UltraSheer	9

BSN Medical Inc.

4455 Autoroute Laval West, Suite 255
Laval (QC) Canada H7P 4W6

www.jobstcanada.ca
www.bsnmedical.ca
Tel 1-877-978-5526
Fax 1-877-978-9703

PE161027RET R2

